


WAY-359473

Chemical Properties

CAS No. : 314248-20-9
 Formula: C₁₄H₁₇N₃S
 Molecular Weight: 259.37
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-359473 is an active compound.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8555 mL	19.2775 mL	38.555 mL
5 mM	0.7711 mL	3.8555 mL	7.711 mL
10 mM	0.3855 mL	1.9277 mL	3.8555 mL
50 mM	0.0771 mL	0.3855 mL	0.7711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481