

WAY-358981

Chemical Properties


CAS No. : 307322-32-3

Formula: C₁₄H₁₂N₄O

Molecular Weight: 252.27

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-358981 is a high-purity biochemical reagent suitable for biochemical experiments and drug syntheses in research.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 20 mg/mL (79.28 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.964 mL	19.820 mL	39.6401 mL
5 mM	0.7928 mL	3.964 mL	7.928 mL
10 mM	0.3964 mL	1.982 mL	3.964 mL
50 mM	0.0793 mL	0.3964 mL	0.7928 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Eglen R M, Bley K, Bonhaus D W, et al. RS 23597-190: a potent and selective 5-HT₄ receptor antagonist[J]. British journal of pharmacology, 1993, 110(1): 119-126.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481