


WAY-313318

Chemical Properties

CAS No. : 351993-88-9
 Formula: C₁₃H₁₀Cl₂N₄
 Molecular Weight: 293.15
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-313318 is an active compound.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4112 mL	17.0561 mL	34.1122 mL
5 mM	0.6822 mL	3.4112 mL	6.8224 mL
10 mM	0.3411 mL	1.7056 mL	3.4112 mL
50 mM	0.0682 mL	0.3411 mL	0.6822 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481