

Vonsetamig

Chemical Properties

CAS No. : 2408319-33-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Vonsetamig, a humanized immunoglobulin G4-kappa monoclonal antibody, targets both TNFRSF17 and CD3E. It serves as an antineoplastic agent.
Targets(IC50)	Integrin,TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481