

## Vitronectin

## Chemical Properties

CAS No. : 83380-82-9

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Vitronectin, a multifunctional glycoprotein found in blood and the extracellular matrix, interacts with glycosaminoglycans, collagen, plasminogen, and the urokinase-receptor, in addition to stabilizing the inhibitory conformation of plasminogen activation inhibitor-1. It is utilized in research related to wound healing and tumor progression [1].
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481