

UR-MB-355

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	UR-MB-355, a GPR3 ligand, demonstrates submicromolar affinities for its primary target GPR3 and the closely related receptors GPR6 and GPR12 [1].
Targets(IC50)	Others
In vitro	UR-MB-355 is a ligand for GPR3, which also binds with the closely related receptors GPR6 and GPR12, exhibiting similar submicromolar affinity.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481