

## Tyrosinase-IN-15

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Tyrosinase-IN-15 (Compound 39), with an IC50 of 7.12 $\mu$ M and a Ki of 11.8 $\mu$ M, functions as a tyrosinase inhibitor [1].
Targets(IC50)	Tyrosinase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481