

Tuparstobart

Chemical Properties

CAS No. : 2649466-03-3

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tuparstobart (Incagn-02385) is an IgG1k monoclonal antibody targeting the immune checkpoint receptor protein LAG-3, predominantly expressed on activated T cells, NK cells, B cells, and plasmacytoid dendritic cells [1].
Targets(IC50)	LAG-3

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481