

## Torapsel

## Chemical Properties

CAS No. : 204658-47-9

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Torapsel (42-89-glycoprotein; WAY164339), a fusion protein comprising immunoglobulin, is utilized for the investigation of ischemia-reperfusion injury prophylaxis.
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481