

## Thyrotropin

## Chemical Properties

CAS No. : 9002-71-5

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Thyrotropin (TSH), a thyroid-stimulating hormone, is synthesized by thyrotrope cells in the anterior pituitary gland and modulates the endocrine activity of the thyroid [1] [2].
Targets(IC50)	Apoptosis,Others,PKC,TSH Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481