

Talizumab

Chemical Properties

CAS No. : 380610-22-0

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Talizumab (TNX 901) is a humanized IgG1 monoclonal antibody that specifically targets circulating immunoglobulin E (IgE). It binds to free IgE not already bound by the high-affinity IgE receptors (FcεRI) present on mast cells and basophils, which leads to the marked downregulation of these receptors and additionally inhibits allergen-specific T-cell activation by interfering with antigen processing and presentation activities.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481