

Sucantomotide

Chemical Properties

CAS No. : 2378651-62-6


Formula: C₅₄H₈₄N₁₀O₁₂S

Molecular Weight: 1097.37

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Sucantomotide is an immunological agent employed for active immunization with antineoplastic properties [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9113 mL	4.5563 mL	9.1127 mL
5 mM	0.1823 mL	0.9113 mL	1.8225 mL
10 mM	0.0911 mL	0.4556 mL	0.9113 mL
50 mM	0.0182 mL	0.0911 mL	0.1823 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481