


SPPHPSPAFSPAFDNLYYWDQ

Chemical Properties

CAS No. : 1443436-80-3
 Formula: C116H150N26O33
 Molecular Weight: 2436.59
 Storage: Keep away from moisture
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SPPHPSPAFSPAFDNLYYWDQ is a multi-epitope class II rat HER2/neu peptide used in the investigation of class II HER2-DC1 vaccine studies [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4104 mL	2.052 mL	4.1041 mL
5 mM	0.0821 mL	0.4104 mL	0.8208 mL
10 mM	0.041 mL	0.2052 mL	0.4104 mL
50 mM	0.0082 mL	0.041 mL	0.0821 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481