

SMP-33693

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SMP-33693 is an Antibody-Drug Conjugate (ADC) known for its low payload shedding rate and demonstrated in vivo tumor inhibition in ovarian, gastric, and breast cancers [1].
Targets(IC50)	Drug-Linker Conjugates for ADC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481