

SAR441255

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	SAR441255 is a potent unimolecular peptide and GLP-1/GIP/GCG receptor triagonist, demonstrating balanced activation across all three target receptors, with positive acute glucoregulatory effects in diabetic obese monkeys [1].
Targets(IC50)	Glucagon Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481