

## Relaxin H3 (human)

## Chemical Properties

CAS No. : 1158181-62-4

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Relaxin H3 (human), a relaxin peptide, exhibits antifibrotic effects by activating RXFP1 [1].
Targets(IC50)	RXFP receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481