

Puxitatug

Chemical Properties

CAS No. : 2760250-47-1

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Puxitatug (INT-016) is a humanised monoclonal antibody targeting VTCN1/B7H4, commonly used as the targeting component of the antibody-drug conjugate (ADC) Puxitatug samrotecan, delivering the drug specifically to cancer cells expressing B7-H4.
Targets(IC50)	Immunology/Inflammation related,B7

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481