

## pTH-Related Protein (1-37) (human, mouse, rat)

### Chemical Properties

CAS No. : 206010-80-2

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	PTH-Related Protein (1-37) (human, mouse, rat) is a peptide analog of parathyroid-hormone-related protein (pTHrP), which plays a role in the physiological regulation of bone formation [1].
Targets(IC50)	Thyroid hormone receptor(THR)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481