

pTH (2-38) (human)

Chemical Properties

CAS No. : 154765-04-5

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

PTH (2-38) (human) plays a role in bone anabolism and can elevate serum levels of INS-PTH (complete N-terminal specific PTH) [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481