

## Proteasome-IN-5

## Chemical Properties


CAS No. : 1613134-34-1

Formula: C<sub>20</sub>H<sub>30</sub>BN<sub>5</sub>O<sub>7</sub>

Molecular Weight: 463.29

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Proteasome-IN-5, also known as compound 5, acts as a proteasome inhibitor [1].
Targets(IC50)	Proteasome

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1585 mL	10.7924 mL	21.5848 mL
5 mM	0.4317 mL	2.1585 mL	4.317 mL
10 mM	0.2158 mL	1.0792 mL	2.1585 mL
50 mM	0.0432 mL	0.2158 mL	0.4317 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481