


PROTAC IRAK4 ligand-3

Chemical Properties

CAS No. :	2434848-46-9
Formula:	C ₂₃ H ₂₅ F ₂ N ₇ O ₃
Molecular Weight:	485.49
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTAC IRAK4 Ligand-3 serves as a ligand for PROTAC IRAK4 Degrader-7 and is utilized in cancer research [1].
Targets(IC50)	Ligands for Target Protein for PROTAC,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0598 mL	10.2989 mL	20.5977 mL
5 mM	0.412 mL	2.0598 mL	4.1195 mL
10 mM	0.206 mL	1.0299 mL	2.0598 mL
50 mM	0.0412 mL	0.206 mL	0.412 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481