


Procyanidin A5'

Chemical Properties

CAS No. :	111466-30-9
Formula:	C30H24O12
Molecular Weight:	576.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Procyanidin A5 is a type A dimerized proanthocyanidin [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7346 mL	8.673 mL	17.3461 mL
5 mM	0.3469 mL	1.7346 mL	3.4692 mL
10 mM	0.1735 mL	0.8673 mL	1.7346 mL
50 mM	0.0347 mL	0.1735 mL	0.3469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481