


Praeroside IV

Chemical Properties

CAS No. : 117233-34-8
 Formula: C₂₀H₂₄O₉
 Molecular Weight: 408.4
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Praeroside IV, a compound extractable from the roots of <i>Angelica furcijuga</i> , is utilized in the research of D-GalN/LPS-induced liver injury [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4486 mL	12.2429 mL	24.4858 mL
5 mM	0.4897 mL	2.4486 mL	4.8972 mL
10 mM	0.2449 mL	1.2243 mL	2.4486 mL
50 mM	0.049 mL	0.2449 mL	0.4897 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481