

PITCOIN4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PITCOIN4 is a Class II Alpha PI3K-C2 α inhibitor with high selectivity, demonstrating nanomolar inhibition of PI3K-C2 α and exhibiting over 100-fold selectivity against a broad kinase panel [1].
Targets(IC50)	PI3K

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481