

Piscerygenin

Chemical Properties


CAS No. : 67978-85-2

Formula: C₁₆H₁₂O₇

Molecular Weight: 316.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Piscerygenin (compound 4), a naturally derived isoflavonoid, can be isolated from the root bark of <i>Piscidia erythrina</i> [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.162 mL	15.8098 mL	31.6196 mL
5 mM	0.6324 mL	3.162 mL	6.3239 mL
10 mM	0.3162 mL	1.581 mL	3.162 mL
50 mM	0.0632 mL	0.3162 mL	0.6324 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481