

Pidilizumab

Chemical Properties

CAS No. : 1036730-42-3

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pidilizumab (CT-011) is a humanized monoclonal antibody that modulates anti-tumor immune responses by targeting the PD-1 receptor and the Notch ligand DLL1. It enhances the activity of T cells and natural killer (NK) cells, demonstrating significant clinical potential in the study of hematologic malignancies, particularly follicular lymphoma.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481