

Pegsebrenatide

Chemical Properties

CAS No. : 2243292-26-2

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pegsebrenatide (NLY01), a long-acting GLP-1R agonist, exhibits an extended half-life and favorable blood-brain barrier penetration. It impedes A1 astrocyte transformation, mitigates dopaminergic cell death, and ameliorates motor symptoms in mouse models of Parkinson's disease (PD) [1].
Targets(IC50)	Glucagon Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481