

Pegloxenatide

Chemical Properties

CAS No. : 2420483-82-3

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pegloxenatide, a glucagon-like peptide-1 receptor (GLP-1 RA) agonist, is utilized in type 2 diabetes research [1].
Targets(IC50)	Apoptosis,Akt,AMPK,Autophagy,Glucagon Receptor,Interleukin,PI3K,ROS,Sirtuin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481