

Patecibart

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Patecibart, a humanized immunoglobulin G4-kappa monoclonal antibody, acts as an antagonist to the endothelin receptor A (EDNRA).
Targets(IC50)	Endothelin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481