

PAP 248-286

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Prostatic Acid Phosphatase (248-286), or PAP (248-286), is a biologically active peptide and a semen-derived enhancer of viral infection (SEVI) factor present in semen. It significantly amplifies HIV infection by promoting enhanced virion attachment to target cells.
Targets(IC50)	HIV Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481