

## Pancreastatin (swine)

## Chemical Properties

CAS No. : 106477-83-2

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Pancreastatin (swine), a 49-residue peptide derived from porcine pancreas [1], significantly suppresses the release of insulin in response to glucose.
Targets(IC50)	Serine Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481