

Palmitoyl tripeptide-8

Chemical Properties

CAS No. : 936544-53-5


Formula: C₃₇H₆₁N₉O₄

Molecular Weight: 695.94

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Palmitoyl tripeptide-8, a bioactive peptide known for its anti-allergen properties, is reported to be used as a cosmetic ingredient [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4369 mL	7.1845 mL	14.3691 mL
5 mM	0.2874 mL	1.4369 mL	2.8738 mL
10 mM	0.1437 mL	0.7185 mL	1.4369 mL
50 mM	0.0287 mL	0.1437 mL	0.2874 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481