

Palmitoyl tetrapeptide-20

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Palmitoyl Tetrapeptide-20 (PTP20) is a bioactive peptide known to preserve the function of follicular melanocytes and is reported to be used as an ingredient in cosmetics [1].
Targets(IC50)	Melanocortin Receptor,Sirtuin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481