


Palmitoyl dipeptide-5 diaminohydroxybutyrate

Chemical Properties

CAS No. :	794590-33-3
Formula:	C ₃₁ H ₆₁ N ₅ O ₅
Molecular Weight:	583.85
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Palmitoyl Dipeptide-5 Diaminohydroxybutyrate is a bioactive peptide recognized for its anti-aging properties and is commonly utilized as an ingredient in cosmetics [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7128 mL	8.5638 mL	17.1277 mL
5 mM	0.3426 mL	1.7128 mL	3.4255 mL
10 mM	0.1713 mL	0.8564 mL	1.7128 mL
50 mM	0.0343 mL	0.1713 mL	0.3426 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481