


Ophiopogonone C

Chemical Properties

CAS No. :	477336-77-9
Formula:	C ₁₉ H ₁₄ O ₇
Molecular Weight:	354.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ophiopogonone C (compound 3), an isoflavone, has been isolated from the ethanol extract of the tubers of <i>Ophiopogon japonicus</i> [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8224 mL	14.1119 mL	28.2239 mL
5 mM	0.5645 mL	2.8224 mL	5.6448 mL
10 mM	0.2822 mL	1.4112 mL	2.8224 mL
50 mM	0.0564 mL	0.2822 mL	0.5645 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481