

## Oligopeptide-74

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Oligopeptide-74, a bioactive peptide known for promoting hair growth, is utilized as a cosmetic ingredient [1].
-------------	---

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481