

Oligopeptide-68

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Oligopeptide-68, a bioactive peptide with anti-aging properties, is commonly utilized as an ingredient in cosmetics [1]. |
| Targets(IC50) | Amino Acids and Derivatives |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481