


Oligopeptide-24

Chemical Properties

CAS No. :	1018685-16-9
Formula:	C50H78N16O19S2
Molecular Weight:	1271.38
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Oligopeptide-24 (CG-EDP3) is a bioactive peptide recognized for its skin repair properties, frequently utilized as an ingredient in cosmetics [1].
-------------	----------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7865 mL	3.9327 mL	7.8655 mL
5 mM	0.1573 mL	0.7865 mL	1.5731 mL
10 mM	0.0787 mL	0.3933 mL	0.7865 mL
50 mM	0.0157 mL	0.0787 mL	0.1573 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481