


## Oligopeptide-10

## Chemical Properties

CAS No. :	466691-40-7
Formula:	C78H140N20O15
Molecular Weight:	1598.07
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Oligopeptide-10, a bioactive peptide known for its antimicrobial properties, is utilized as an ingredient in cosmetics [1].
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6258 mL	3.1288 mL	6.2575 mL
5 mM	0.1252 mL	0.6258 mL	1.2515 mL
10 mM	0.0626 mL	0.3129 mL	0.6258 mL
50 mM	0.0125 mL	0.0626 mL	0.1252 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481