


OAT1/3-IN-1

Chemical Properties

CAS No. :	38393-73-6
Formula:	C ₁₆ H ₁₂ O ₆
Molecular Weight:	300.26
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	OAT1/3-IN-1 (compound 7), a dual inhibitor of OAT1 and OAT3, mitigates Cys-Hg toxicity in HEK-OAT1 cells at 10 μ M and offers potential kidney protection. This compound is used in the investigation of mercury-induced renal damage [1].
Targets(IC50)	OAT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3304 mL	16.6522 mL	33.3045 mL
5 mM	0.6661 mL	3.3304 mL	6.6609 mL
10 mM	0.333 mL	1.6652 mL	3.3304 mL
50 mM	0.0666 mL	0.333 mL	0.6661 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481