

NLRP3-IN-23

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | NLRP3-IN-23 (Compound 15C) is an inhibitor that significantly impedes heme-mediated activation of the NLRP3 inflammasome at a concentration of 0.1 μ M [1]. |
| Targets(IC50) | NOD-like Receptor (NLR) |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481