

## Neuromedin S(rat) TFA

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Neuromedin S(rat) TFA, a 34-amino acid peptide derivative of rat Neuromedin S, is a neuropeptide sourced from the rat brain that functions as a ligand for the G protein-coupled receptor FM4/TGR-1 [1].
-------------	--

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481