

Neuromedin S (human)

Chemical Properties

CAS No. : 1138204-27-9

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Neuromedin S (human), a 33-amino acid neuropeptide, serves as an endogenous ligand for the orphan G-protein coupled receptor (GPCR) FM-4/TGR-1, and acts on the neuromedin U (NMU) receptor 2 (NMUR2) to regulate body weight homeostasis [1].
Targets(IC50)	NMUR, GPCR19

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481