


MPS-Gαi3

Chemical Properties

CAS No. : 439223-21-9
Formula: C120H205N29O31S
Molecular Weight: 2582.15
Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	MPS-Gαi3, a cell-penetrating peptide, facilitates drug delivery [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3873 mL	1.9364 mL	3.8727 mL
5 mM	0.0775 mL	0.3873 mL	0.7745 mL
10 mM	0.0387 mL	0.1936 mL	0.3873 mL
50 mM	0.0077 mL	0.0387 mL	0.0775 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481