


MBP (146-170)

Chemical Properties

CAS No. :	604803-17-0
Formula:	C112H194N40O34S
Molecular Weight:	2677.05
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Myelin Basic Protein (MBP) fragment 146-170 is utilized in the study of Multiple Sclerosis (MS) [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3735 mL	1.8677 mL	3.7355 mL
5 mM	0.0747 mL	0.3735 mL	0.7471 mL
10 mM	0.0374 mL	0.1868 mL	0.3735 mL
50 mM	0.0075 mL	0.0374 mL	0.0747 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481