

Maridebart

Chemical Properties

CAS No. : 2761478-99-1

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Maridebart is a fully humanized monoclonal antibody targeting gastric inhibitory polypeptide receptor (GIPR), and it is also the antibody part of Maridebart cafraglutide, which can be used to construct antibody-drug conjugate (ADC). Maridebart is mainly used to study the mechanism of metabolic diseases such as obesity and type 2 diabetes.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481