

MAQAAEYR

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	MAQAAEYR, an antioxidant peptide derived from the pancreatin hydrolysate of C-phycocyanin, exhibits potent scavenging activities against DPPH and ABTS radicals. It also confers protection to zebrafish larvae against oxidative damage induced by H <sub>2</sub> O <sub>2</sub> [1].
Targets(IC50)	Reactive Oxygen Species

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481