

Maculatin 1.1 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Maculatin 1.1 TFA, an antimicrobial peptide, exhibits a minimum inhibitory concentration (MIC) of 7 μ M against Staphylococcus aureus and induces bacterial mortality by perforating the organism's membrane [1].
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481