


M351-110

Chemical Properties

CAS No. : 1215797-40-2
 Formula: C₁₅H₁₂ClFN₂O₂S₃
 Molecular Weight: 402.91
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	M351-110 is an agonist of the V-domain Ig suppressor of T cell activation (VISTA) designed for oncological studies [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4819 mL	12.4097 mL	24.8194 mL
5 mM	0.4964 mL	2.4819 mL	4.9639 mL
10 mM	0.2482 mL	1.241 mL	2.4819 mL
50 mM	0.0496 mL	0.2482 mL	0.4964 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481