


## L-Tryptophanyl-glycine

## Chemical Properties

CAS No. :	7360-09-0
Formula:	C <sub>13</sub> H <sub>15</sub> N <sub>3</sub> O <sub>3</sub>
Molecular Weight:	261.28
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	L-Tryptophanyl-glycine (Trp-gly), a dipeptide, inhibits Ceronapril uptake in Caco-2 cells with an inhibition constant (K <sub>i</sub> ) of 7.19 mM [1].
Targets(IC <sub>50</sub> )	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8273 mL	19.1366 mL	38.2731 mL
5 mM	0.7655 mL	3.8273 mL	7.6546 mL
10 mM	0.3827 mL	1.9137 mL	3.8273 mL
50 mM	0.0765 mL	0.3827 mL	0.7655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481